



TECHNICAL LITERATURE COMPACT SERIES BUILT IN VERSIONS

The **Compact** Built In series of batteries provides a full days power for most users. Batteries are available for the latest Otto Bock™ adult powered components, VASi components and all RSLSTEEPER powered systems. In addition, a 5Volt version can be supplied for use with Otto Bock System 2000 series hands. Built in versions in both voltages are available where use is permitted.

At a glance

- ❑ Advanced safety technology prevents excessive demands on the battery.
- ❑ Compatible with Otto Bock DMC and other hands, wrists and Greifers. This battery can also
- ❑ be used with all RSLSTEEPER hands including the Scamp series.
- ❑ 5-Volt version available for Otto Bock system 2000 series hands.
- ❑ “All Day” capacity (500MAhr).
- ❑ 80% of capacity is restored in the first hour of charging.
- ❑ Low self discharge in storage.

IMPORTANT NOTES: -

1)

! Warnings! :

- **Not to be charged at greater than 1C**
- **Not to be charged when connected to patient**
- **No user serviceable parts inside**
- **When the battery is to be first used from store we recommend that it be charged prior to issue for a minimum of 19 hours.**
- **We recommend that the battery is charged each night so that all day capacity is available**
- **Do not charge for more than 24 hours**
- **Warning: Do not use with Nicel Cadmium or Nickel Metal Hydride battery charges. Use only the RSLSTEEPER Charger supplied for these Li –Ion Batteries 1) Lithium-ion cells and battery packs may get hot, explode or ignite and cause serious injury if exposed to abuse conditions. Be sure to follow the safety warnings listed below:**
 - Do not place the battery in fire or heat the battery.
 - Do not use the battery other than directed by the Prosthetist.

- **Do not connect the positive terminal and negative terminal of the battery to each other with any metal object (such as wire).**
- **Do not pierce the battery with nails, strike the battery with a hammer, step on the battery or otherwise subject it to strong impacts or shocks.**
- **Do not expose battery to water or salt water, or allow the battery to get wet.**

2) Do not disassemble or modify the battery. The battery contains safety and protection devices, which, if damaged, may cause the battery to generate heat, explode or ignite.

3) If the battery case becomes damaged stop using and place in a plastic bag to stop the battery getting wet and return to prosthetist.

4) Do not place the battery in or near fire, on stoves or other high temperature locations. Do not place the battery in direct sunlight, or use or store the battery inside cars in hot weather. Doing so may cause the battery to generate heat, explode or ignite. Using the battery in this manner may also result in a loss of performance and a shortened life expectancy.

5) In the unlikely event of a malfunction or unusual operation stop using and consult the prosthetist

! CAUTION!

6) If the device is to be used by small children, the prosthetist / caregiver should explain the contents of this document to the children and provide adequate supervision to ensure the device is being used appropriately.

7) Immediately discontinue use of the battery if, while using, charging or storing the battery, the battery emits an unusual smell, feels hot, changes color or shape, or appears abnormal in any other way. Contact RSLSTEEPER if any of these problems are observed.

8) Do not place the battery in microwave ovens, high-pressure containers or on induction cookware.

9) In the event the battery is damaged which results in the cell leaking and the fluid gets into one's eye, do not rub the eye. Rinse wells with water and immediately seek medical care. If left untreated, the battery fluid could cause damage to the eye.

For further information contact:

RSLSTEEPER, Products Division Unit 7, Hunslet Trading Estate Severn Road Leeds West
YorkshireLS10 1BL United Kingdom
Tel: +44 (0) 870 240 4133 / Fax: +44 (0) 0870 230 4779 Email: sales@rehab.co.uk
www.rslsteeper.com

The design and manufacture of RSL STEEPER equipment and components are subject to a policy of continuous reappraisal. The Company therefore reserves the right to introduce changes and withdraw products without notice. © COPYRIGHT 2003,

RSL STEEPER. Technical Literature L23348 issue 1 April 2003 