

# ***RSLSteeper ARGO 60***

*weight limit 60 Kilos (132 pounds)*



*ARGO 60 Assembled*

## **AR300 Basic Intermediate Advanced Reciprocating Gait Orthosis Kit**

- *for small adults and teen-age patients weighing up to 60 kilos (132 pounds)*
- *shown assembled*



- AR311** Back Tube, X-Small (brown) use Drive Cable AR101
- AR312** Back Tube, Small (blue) use Drive Cable AR101
- AR313** Back Tube, Medium-Small (purple) use Drive Cable AR102
- AR314** Back Tube, Medium (red) use Drive Cable AR102
- AR315** Back Tube, Medium-Large (yellow) use Drive Cable AR102
- AR316** Back Tube, Large (black) use Drive Cable AR103
- AR317** Back Tube, X-Large (green) use Drive Cable AR103



- AR101** Drive Cable Kit, Small (375 mm, 14.8")
- AR102** Drive Cable Kit, Medium-Small (415 mm, 16.3")
- AR105** Drive Cable Kit, Medium (455 mm, 17.9")
- AR103** Drive Cable Kit, Medium-Large (495 mm, 19.5")
- AR104** Drive Cable Kit, Large (555 mm, 21.9")



- AR321** Knee-Lock Cables, pair, 165 mm, brown
- AR322** Knee-Lock Cables, pair, 195 mm, blue
- AR323** Knee-Lock Cables, pair, 225 mm, purple
- AR324** Knee-Lock Cables, pair, 255 mm, red
- AR325** Knee-Lock Cables, pair, 285 mm, yellow
- AR326** Knee-Lock Cables, pair, 315 mm, black
- AR327** Knee-Lock Cables, pair, 345 mm, green
- AR328** Knee-Lock Cables, pair, 375 mm, white

## *RSLSteeper ARGO 60*



**AR300A-KO** Abduction Hip Joint Kit, No Knee Assist  
**AR300A-K2** Abduction Hip Joint Kit, 200N Knee Assist  
**AR300A-K3** Abduction Hip Joint Kit, 300N Knee Assist  
**AR300A-K4** Abduction Hip Joint Kit, 400N Knee Assist



**AR300F-KO** Fixed Hip Joint Kit, No Knee Assist  
**AR300F-K2** Fixed Hip Joint Kit, 200N Knee Assist  
**AR300F-K3** Fixed Hip Joint Kit, 300N Knee Assist  
**AR300F-K4** Fixed Hip Joint Kit, 400N Knee Assist



Extension-Assist Knee Joint options:

**K2** 91 pounds (200N)  
**K3** 136 pounds (300N)  
**K4** 182 pounds (400N)



**AR161** Universal Foot Plates, pair, women's - size 8

Use with Becker modular double-action ankle joints (2825-A)