

RSLSteeper ARGO 90

weight limit 90 Kilos (200 pounds)



ARGO 90 Assembled

AR100 - Basic Adult Advanced Reciprocating Gait Orthosis Kit ARGO90

- *for adult patients weighing up to 90 kilos (200 pounds)*
- *shown assembled*



- AR118** Back Tube, X-Small (white) use Drive Cable AR101
- AR106** Back Tube, Small (blue) use Drive Cable AR101
- AR119** Back Tube, Medium-Small (purple) use Drive Cable AR102
- AR107** Back Tube, Medium (red) use Drive Cable AR102
- AR108** Back Tube, Medium-Large (yellow) use Drive Cable AR102
- AR109** Back Tube, Large (black) use Drive Cable AR103
- AR116** Back Tube, X-Large (green) use Drive Cable AR103




- AR101** Drive Cable Kit, Small (375 mm, 14.8")
- AR102** Drive Cable Kit, Medium-Small (415 mm, 16.3")
- AR105** Drive Cable Kit, Medium (455 mm, 17.9")
- AR103** Drive Cable Kit, Medium-Large (495 mm, 19.5")
- AR104** Drive Cable Kit, Large (555 mm, 21.9")



- AR117** Knee-Lock Cables, pair, 180 mm, white
- AR115** Knee-Lock Cables, pair, 210 mm, yellow-green
- AR110** Knee-Lock Cables, pair, 240 mm, yellow
- AR111** Knee-Lock Cables, pair, 270 mm, blue (shown)
- AR112** Knee-Lock Cables, pair, 300 mm, black
- AR113** Knee-Lock Cables, pair, 330 mm, red

RSLSteeper ARGO 90



AR160 Universal Foot Plates, pair, men's - size 10
AR161 Universal Foot Plates, pair, women's - size 8

Use with Becker Orthopedic modular double-action ankle joints (2825-A)



AR120 Growth Kit, (new L/R Upper Steels)

Requires custom fitting by LTI



AR123 Upper Steel, Extension-Assist Knee, left (shown)
AR124 Upper Steel, Extension-Assist Knee, right

Heavy-Duty versions:

AR123H Heavy-Duty Upper Steel, Extension-Assist Knee, left
AR124H Heavy-Duty Upper Steel, Extension-Assist Knee, right