

Battery Mounts & Adapters

Battery Mounts



BM92 Holder, Compact Battery, Standard Depth, Bock Connector
BM94 Holder, Compact Battery, Standard Depth, LTI Connector
For use with Li-ion Compact batteries BP92 & BP93. Requires opening measuring 1.1"x2.56" (28x65mm). Holder penetrates 0.75" (19mm) below surface. Installed height to center of battery is 0.87" (22mm) above bottom of Holder. Requires one of the MC30 series Battery Adapters.



BM93 Holder, Compact Battery, Shallow Depth, Bock Connector
BM95 Holder, Compact Battery, Shallow Depth, LTI Connector
For use with Li-ion Compact batteries BP92 & BP93. Requires opening measuring 1.1"x2.56" (28x65mm). Holder penetrates 0.57" (14.5mm) below surface. Installed height to center of battery is 0.87" (22mm) above bottom of Holder. Requires one of the MC30 series Battery Adapters.



Adapters for use with Recharge/ON-Off Modules

MC30 Battery Adapter, Utah ProControl Plug



MC31 Battery Adapter, Bock QD Plug



MC34 Battery Adapter, Steeper Plug



Battery Mounts & Adapters

MC37 Battery Adapter, one Bock 13E125, VASI Plug



VASI Plug



MC38 Battery Adapter, two Bock 13E125,



MC67 Battery Adapter, one Bock 13E125, Steeper Plug



MC68 Battery Adapter, two Bock 13E125, Steeper Plug

