

## ***Battery Chargers***

### ***Battery Chargers for LTI Built-in Batteries:***

The battery chargers below are to be used with LTI Built-in Batteries. They connect to the battery through one of the Recharge On-Off Modules. The barrel connectors on these Recharge Modules are sized to mate with the appropriate battery charger to assure that the charger's output is appropriate for the battery. Always use the Recharge Module specified for the particular battery.



**BC96** Fast Charger for LTI Large Capacity Lithium Built-in Batteries  
*Use with 750-900 mAhr, 7.4 V batteries – BP88, BP98, BP102*  
*Barrel connector compatible with BP011 & BP018 Recharge/On/Off Switch Modules*



**BC100** Fast Charger for LTI Small Capacity Lithium Built-in Batteries  
*Use with 250-450 mAhr, 7.4 V batteries – BP101 & BP103*  
*Barrel connector compatible with BP013 & BP023 Recharge/On/Off Switch Modules*



**BE366** Fast Charger for LTI Boston Digital Arm System  
*Use with BE360 & BE360+ Lithium-polymer 11 V batteries*

### ***Battery Chargers for Removable Batteries:***

Battery chargers for removable batteries are table-top units that charge two batteries at a time. Indicator lights tell the state of charge of each battery. It is recommended that both batteries be charged overnight to assure that the user has adequate battery capacity to operate their prosthesis for a full day. Batteries are not to be stored in the charger for long periods of time.



**BC72** Fast Charger for BP72 Removable Batteries  
*Use only with the BP72 Removable Lithium-polymer Battery*

## ***Battery Chargers***



### **BC90** Fast Charger for BP90 Removable Batteries

*Use only with the BP90 Removable Lithium-polymer Battery*

*Do not use older Ni-Cad chargers with the BP90 Battery!*

*Do not use Bock 757B8 Ni-CAD Batteries with this charger!*