

# SKINERGY Leg Covers

## SKINERGY™ Production Leg Covers:

SKINERGY™ protective production leg covers are manufactured using Tufsil 2, the latest innovation in silicone polymer technology. SKINERGY provides an extremely high level of cosmesis for the below knee prostheses, incorporating a unique seamless finish. Its improved tear resistance and 25% elasticity accommodates virtually all sockets and protects the prosthetic components of the leg. And each cover stretches to cover two foot sizes. SKINERGY is also stain resistant for longer life and can easily be cleaned with soap and water.

SKINERGY is available in both split-toe and closed-toe versions making it suitable for all manufacturers' prosthetic feet and with two calf sizes, it can accommodate all users. Both styles are available in 18 colors with nails are hand painted for a "semi-custom" finish. Colors are selected from a color swatch sample. SKINERGY can be rolled down to access the underlying prosthetic hardware for maintenance and adjustments. To access the foot bolt for maintenance purposes, a 25 mm punch tool is available. And if the leg cover is damaged, a patch kit is offered consisting of color-matched patches and adhesive. SKINERGY leg covers are available in four foot sizes and two calf sizes. The foot size is designated by the number (i.e. 22) which relates to the foot size in cm (22 cm or 220 mm). Next is the side which is specified by either a "L" or "R". Finally, each foot size has two calf sizes "small" or "large" which is designated by the trailing letter (S or L).

## Split-Toe Skinergy:



**TLC2XXX SKINERGY Protective Leg Cover, Split Toe**  
For use with all split-toe prosthetic feet

Part Number	Foot Size	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*	
		minimum	maximum	minimum	maximum	minimum	maximum	Mid Calf	Cover Top
	mm	cm	cm	cm	cm	cm	cm	cm	cm
<b>TLC22LS</b>	220/230	20	25	28	35	26	32.5	32	55
<b>TLC22RS</b>	220/230	20	25	28	35	26	32.5	32	55
<b>TLC22LL</b>	220/230	20	25	32	40	30	37.5	32	55
<b>TLC22RL</b>	220/230	20	25	32	40	30	37.5	32	55
<b>TLC24LS</b>	240/250	22	27.5	31	39	29	36	34.5	57
<b>TLC24RS</b>	240/250	22	27.5	31	39	29	36	34.5	57
<b>TLC24LL</b>	240/250	22	27.5	35	44	33	41	34.5	57
<b>TLC24RL</b>	240/250	22	27.5	35	44	33	41	34.5	57
<b>TLC26LS</b>	260/270	22	28	32	40	30	37.5	37	61
<b>TLC26RS</b>	260/270	22	28	32	40	30	37.5	37	61
<b>TLC26LL</b>	260/270	22	28	37	46	35	44	37	61
<b>TLC26RL</b>	260/270	22	28	37	46	35	44	37	61
<b>TLC28LS</b>	280/300	24	30	34	42.5	32	40	39.5	66
<b>TLC28RS</b>	280/300	24	30	34	42.5	32	40	39.5	66
<b>TLC28LL</b>	280/300	24	30	39	49	37	46	39.5	66
<b>TLC28RL</b>	280/300	24	30	39	49	37	46	39.5	66

# SKINERGY Leg Covers

## Closed-Toe Skinergy:



**TNS2XXX SKINERGY Protective Leg Cover, Closed Toe**  
*For use with all closed-toe prosthetic feet*

Part	Foot	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*	
Number	Size	minimum	maximum	minimum	maximum	minimum	maximum	Mid Calf	Cover Top
	mm	cm	cm	cm	cm	cm	cm	cm	cm
<b>TNS22LS</b>	220/230	20	25	28	35	26	32.5	32	55
<b>TNS22RS</b>	220/230	20	25	28	35	26	32.5	32	55
<b>TNS22LL</b>	220/230	20	25	32	40	30	37.5	32	55
<b>TNS22RL</b>	220/230	20	25	32	40	30	37.5	32	55
<b>TNS24LS</b>	240/250	22	27.5	31	39	29	36	34.5	57
<b>TNS24RS</b>	240/250	22	27.5	31	39	29	36	34.5	57
<b>TNS24LL</b>	240/250	22	27.5	35	44	33	41	34.5	57
<b>TNS24RL</b>	240/250	22	27.5	35	44	33	41	34.5	57
<b>TNS26LS</b>	260/270	22	28	32	40	30	37.5	37	61
<b>TNS26RS</b>	260/270	22	28	32	40	30	37.5	37	61
<b>TNS26LL</b>	260/270	22	28	37	46	35	44	37	61
<b>TNS26RL</b>	260/270	22	28	37	46	35	44	37	61
<b>TNS28LS</b>	280/300	24	30	34	42.5	32	40	39.5	66
<b>TNS28RS</b>	280/300	24	30	34	42.5	32	40	39.5	66
<b>TNS28LL</b>	280/300	24	30	39	49	37	46	39.5	66
<b>TNS28RL</b>	280/300	24	30	39	49	37	46	39.5	66

Dimensions are in cm, except foot size which is in mm

\* Cover Length – from floor to mid-calf or cover top



**P80 Cosmetic Cover Lubricant, Pint**  
*makes installing leg covers easier, leaves no residue*