

Cosmetic & Protective Covers

Custom Silicone Cosmetic Covers & Restorations



- High-Definition
- Infinite colors for optimal match
- Stain-resistant silicone
- Silicone or acrylic finger/toe nails
- Match user's anatomy
- Optional features include: hair, freckles, veins and tattoos
- 5-6 week delivery

High-Definition Custom Silicone Leg Cover

RSLSteeper's custom silicone cosmetic covers are produced using a unique fabrication method. Unlike most cosmetic covers that are produced with clear silicone and then "painted" on the interior surface, these are created by "blending" color pigment into the raw silicone to create a more uniform appearance. Since the color is present throughout the silicone, it has more "depth" and it is not subject to wear and deterioration. Additionally, the density of the silicone can be varied to provide more durability in high-wear areas and softer surfaces at feathered edges where the cosmesis contacts the user's skin. The possibilities are limitless!

High-Definition Silicone covers are the premium cosmesis resulting in a life-like prosthesis. These are custom sculpted with the highest definition to match the user's anatomy. They are created with an infinite number of colors to assure a match with the user's skin tone. Finger nails are crafted out of stiff silicone or acrylic (for nail painting) and use as many colors as is necessary for a match. Optional features include; freckles, special skin coloration, veins, surface hair and even tattoos if required.

*Production cosmetic covers are also available. These products fill the range for appearance, durability and price. **Production PVC** is the standard cosmetic cover used with most prostheses. This is an economical alternative and yet these gloves provide good cosmetic definition and durability. **True Finish™ PVC** has greater detail than standard production PVC and uses a unique multi-color skin pigment for a more natural look. In the silicone line, **Production Silicone** for Otto Bock hands is the least expensive of the silicone product offerings. It has the benefits provided by all silicone covers: greater resiliency and stain resistance. The cosmetic detail approaches that of the standard production PVC gloves.*

Note: 50% deposit required with all custom high-definition silicone cover orders

Cosmetic Restorations - Digits



SC73 Finger/Toe - Premium quality w/acrylic or silicone nail

Cosmetic & Protective Covers

Upper-Extremity Cosmeses



SC82 Partial Hand, 2 fingers - Premium quality
SC83 Partial Hand, 3 fingers - Premium quality



SC85 Full Hand - Premium quality



SC87 Below Elbow - Premium quality



SC89 Above Elbow - Premium quality



SC100 Cosmetic Glove for Otto Bock Hands - Premium quality

Note: 50% deposit required with all custom high-definition silicone cover orders

Cosmetic & Protective Covers

Lower-extremity Cosmeses



SC31 Partial/Complete Foot - Intermediate quality

SC32 Partial/Complete Foot - Premium quality



SC52 Below Knee - Premium quality

SC53 Below Knee for FlexFoot/Total Shock - Premium

Note: 50% deposit required with all custom high-definition silicone cover orders

SKINERGY™ Protective Leg Covers

SKINERGY™ protective leg covers are a new product manufactured using Tufsil 2, the latest innovation in silicone polymer technology. *SKINERGY* provides an extremely high level of cosmesis for the below knee or partial foot prosthesis, incorporating a unique seamless finish. Its improved tear resistance and 25% elasticity accommodates virtually all sockets and protects the prosthetic components of the leg. And each cover stretches to cover two foot sizes. *SKINERGY* is also stain resistant for longer life and can easily be cleaned with soap and water.

SKINERGY is now available in both split-toe and closed-toe versions making it suitable for all manufacturer's prosthetic feet and with two calf sizes, it can accommodate all users. Both styles are available in 12 colors and nails are hand painted for a "semi-custom" finish. Colors are selected from a color swatch sample. *SKINERGY* can be rolled down to access the underlying prosthetic hardware for maintenance and adjustments. To access the foot bolt for maintenance purposes, a 25 mm punch tool is available. And if the leg cover is damaged, a patch kit is offered consisting of color-matched patches and adhesive.

SKINERGY leg covers are available in four foot sizes and two calf sizes. The foot size is designated by the number (i.e. 22) which relates to the foot size in cm (22 cm or 220 mm). Each foot size has two calf sizes "small" or "large" designated by the trailing letter (S or L).

Cosmetic & Protective Covers



LC2XXX SKINERGY Protective Leg Cover, Split Toe
For use with all split-toe prosthetic feet

Part	Foot	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*	
Number	Size mm	minimum cm	maximum cm	minimum cm	maximum cm	minimum cm	maximum cm	Mid Calf cm	Cover Top cm
LC22LS	220/230	20	25	28	35	26	32.5	32	55
LC22RS	220/230	20	25	28	35	26	32.5	32	55
LC22LL	220/230	20	25	32	40	30	37.5	32	55
LC22RL	220/230	20	25	32	40	30	37.5	32	55
LC24LS	240/250	22	27.5	31	39	29	36	34.5	57
LC24RS	240/250	22	27.5	31	39	29	36	34.5	57
LC24LL	240/250	22	27.5	35	44	33	41	34.5	57
LC24RL	240/250	22	27.5	35	44	33	41	34.5	57
LC26LS	260/270	22	28	32	40	30	37.5	37	61
LC26RS	260/270	22	28	32	40	30	37.5	37	61
LC26LL	260/270	22	28	37	46	35	44	37	61
LC26RL	260/270	22	28	37	46	35	44	37	61
LC28LS	280/300	24	30	34	42.5	32	40	39.5	66
LC28RS	280/300	24	30	34	42.5	32	40	39.5	66
LC28LL	280/300	24	30	39	49	37	46	39.5	66
LC28RL	280/300	24	30	39	49	37	46	39.5	66

Dimensions are in cm, except foot size which is in mm

* Cover Length – from floor to mid-calf or cover top



P80 Cosmetic Cover Lubricant, Pint
makes installing leg covers easier, leaves no residue

Cosmetic & Protective Covers



NS2XXX SKINERGY Protective Leg Cover, Closed Toe
For use with all closed-toe prosthetic feet

Part	Foot	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*	
Number	Size mm	minimum cm	maximum cm	minimum cm	maximum cm	minimum cm	maximum cm	Mid Calf cm	Cover Top cm
NS22LS	220/230	20	25	28	35	26	32.5	32	55
NS22RS	220/230	20	25	28	35	26	32.5	32	55
NS22LL	220/230	20	25	32	40	30	37.5	32	55
NS22RL	220/230	20	25	32	40	30	37.5	32	55
NS24LS	240/250	22	27.5	31	39	29	36	34.5	57
NS24RS	240/250	22	27.5	31	39	29	36	34.5	57
NS24LL	240/250	22	27.5	35	44	33	41	34.5	57
NS24RL	240/250	22	27.5	35	44	33	41	34.5	57
NS26LS	260/270	22	28	32	40	30	37.5	37	61
NS26RS	260/270	22	28	32	40	30	37.5	37	61
NS26LL	260/270	22	28	37	46	35	44	37	61
NS26RL	260/270	22	28	37	46	35	44	37	61
NS28LS	280/300	24	30	34	42.5	32	40	39.5	66
NS28RS	280/300	24	30	34	42.5	32	40	39.5	66
NS28LL	280/300	24	30	39	49	37	46	39.5	66
NS28RL	280/300	24	30	39	49	37	46	39.5	66

Dimensions are in cm, except foot size which is in mm

* Cover Length – from floor to mid-calf or cover top



P80 Cosmetic Cover Lubricant, Pint
makes installing leg covers easier, leaves no residue