

Friction Elbows, Elbow Joints & Locking Hinges

Friction Elbows for Children. The RSLSteeper Friction Elbow Units for children provide 135° of flexion and 360° of humeral rotation. Both movements have friction adjustment. The elbows are available in two sizes: a 1½" Elbow (FE150) for children requiring a 2" hand; and a 1¾" Elbow (FE175) for children requiring a 2¼" hand. Both elbow units are supplied with two stainless side steels for forearm attachment. They can be used with passive or electric hands, and have a pass-through for myoelectric hand wires. These elbow joints can also be used as shoulder joints. Manufactured in flesh-pink only.



FE150 1½" Infant's Universal Friction Elbow



FE 175 1¾" Child's Universal Friction Elbow

Friction Elbow Joints



VFE 11 L VASI Small Friction Elbow Joint
VFE 12 R VASI Small Friction Elbow Joint



VFE14 Spare Cover for Friction Elbow Joint

Friction Elbows, Elbow Joints & Locking Hinges



EJ1001 L Outside Heavy-duty Locking Hinge
EJ1002 R Outside Heavy-duty Locking Hinge

*A small (3mm) pull will release the lock. It will re-lock in the next position.
A large (8mm) pull will toggle the mechanism to the fully unlocked state.*



EJ4011 L/R Free-Swing Hinge