

High-Definition Custom Silicone Cosmetic Covers & Restorations



High-Definition Custom Silicone Leg Cover

- **High-Definition**
- **Infinite colors for optimal match**
- **Stain-resistant silicone**
- **Silicone or acrylic finger/toe nails**
- **Match user's anatomy**
- **Optional features include: hair, freckles, veins and tattoos**
- **5-6 week delivery**

RSL Steeper's custom silicone cosmetic covers are produced using a unique fabrication method. Unlike most cosmetic covers that are produced with clear silicone and then "painted" on the interior surface, these are created by "blending" color pigment into the raw silicone to create a more uniform appearance. Since the color is present throughout the silicone, it has more "depth" and it is not subject to wear and deterioration. Additionally, the density of the silicone can be varied to provide more durability in high-wear areas and softer surfaces at feathered edges where the cosmesis contacts the user's skin. The possibilities are limitless!

High-Definition Silicone covers are the premium cosmesis resulting in a life-like prosthesis. These are custom sculpted with the highest definition to match the user's anatomy. They are created with an infinite number of colors to assure a match with the user's skin tone. Finger nails are crafted out of stiff silicone or acrylic (for nail painting) and use as many colors as is necessary for a match. Optional features include; freckles, special skin coloration, veins, surface hair and even tattoos if required.

*Production cosmetic covers are also available. These products fill the range for appearance, durability and price. **Production PVC** is the standard cosmetic cover used with most prostheses. This is an economical alternative and yet these gloves provide good cosmetic definition and durability. **True Finish™ PVC** has greater detail than standard production PVC and uses a unique multi-color skin pigment for a more natural look. In the silicone line, **Production Silicone** for Otto Bock hands is the least expensive of the silicone product offerings. It has the benefits provided by all silicone covers: greater resiliency and stain resistance. The cosmetic detail approaches that of the standard production PVC gloves.*

Note: 50% deposit required with all custom high-definition silicone cover orders

Cosmetic Restorations - Digits



SC73 Finger/Toe - Premium quality w/acrylic or silicone nail

Upper-Extremity Cosmeses



SC82 Partial Hand, 2 fingers - Premium quality
SC83 Partial Hand, 3 fingers - Premium quality



SC85 Full Hand - Premium quality



SC87 Below Elbow - Premium quality



SC89 Above Elbow - Premium quality



SC100 Cosmetic Glove for Otto Bock Hands - Premium quality

Note: 50% deposit required with all custom high-definition silicone cover orders

Lower-extremity Cosmeses



SC31 Partial/Complete Foot - Intermediate quality
SC32 Partial/Complete Foot - Premium quality



SC52 Below Knee - Premium quality
SC53 Below Knee for FlexFoot/Total Shock - Premium

Note: 50% deposit required with all custom high-definition silicone cover orders