

## LTI Data Sheet Select Electric Hands & Gloves

### SELECT HANDS - Quick Disconnect Wrists (QD)

New Hand #	Old* Hand #	at Metacarpal		side	weight	wrist	description	New Select Hand Gloves		Old Multi Control Gloves	
		circ.	width		gms/lb	diameter		Silicone	PVC	LTI - PVC	RSL-PVC
<b>MYA675LQD</b>	H203L-QD	6¾"	2½"	L	320/0.7	1½"/40mm	Childs Hand	TSG340L	TG340L	G51L	CG78
<b>MYA675RQD</b>	H204R-QD	6¾"	2½"	R	320/0.7	1½"/40mm	Childs Hand	TSG344R	TG344R	G52R	CG79
<b>MYA725LQD</b>	H205L-QD	7¼"	2¾"	L	470/1.0	1¾"/45mm	Adolescent	TSG306L	TG306L	G45L	CG80
<b>MYA725RQD</b>	H206R-QD	7¼"	2¾"	R	470/1.0	1¾"/45mm	Adolescent	TSG307R	TG307R	G46R	CG81
<b>MYA725LQD</b>	H205L-QD	7¼"	2¾"	L	470/1.0	1¾"/45mm	Female Adult Hand	TSG304L	TG304L	G45L	CG80
<b>MYA725RQD</b>	H206R-QD	7¼"	2¾"	R	470/1.0	1¾"/45mm	Female Adult Hand	TSG305R	TG305R	G46R	CG81
<b>MYA775LQD</b>	H207L-QD	7½"-7¾"	3"	L	510/1.2	2"/50mm	Female Adult Hand	TSG312L	TG312L	G47L	CG82
<b>MYA775RQD</b>	H208R-QD	7½"-7¾"	3"	R	510/1.2	2"/50mm	Female Adult Hand	TSG313R	TS313R	G48R	CG83
<b>MYA775LQD</b>	H207L-QD	7½"-7¾"	3"	L	510/1.2	2"/50mm	Male Adult Hand	TSG302L	TG302L	G47L	CG82
<b>MYA775RQD</b>	H208R-QD	7½"-7¾"	3"	R	510/1.2	2"/50mm	Male Adult Hand	TSG303R	TG303R	G48R	CG83
<b>MYA825LQD</b>	H209L-QD	7¾"-8¼"	3¼"	L	520/1.4	2⅛"/54mm	Male Adult Hand	TSG300L	TG300L	G49L	CG84
<b>MYA825RQD</b>	H210R-QD	7¾"-8¼"	3¼"	R	520/1.4	2⅛"/54mm	Male Adult Hand	TSG301R	TG301R	G50R	CG85

### SELECT HANDS - North American (½-20) Threaded Stud

New Hand #	Old* Hand #	at Metacarpal		side	weight	wrist	description	New Select Hand Gloves		Old Multi Control Hand Gloves	
		circ.	width		gms/lbs	diameter		Silicone	PVC	LTI - PVC	RSL-PVC
<b>MYA675LNA</b>	H203L-TH	6¾"	2½"	L	303/0.67	1½"/40mm	Childs Hand	TSG340L	TG340L	G51L	CG78
<b>MYA675RNA</b>	H204R-TH	6¾"	2½"	R	303/0.67	1½"/40mm	Childs Hand	TSG344R	TG344R	G52R	CG79
<b>MYA725LNA</b>	H205L-TH	7¼"	2¾"	L	453/1.0	1¾"/45mm	Adolescent	TSG306L	TG306L	G45L	CG80
<b>MYA725RNA</b>	H206R-TH	7¼"	2¾"	R	453/1.0	1¾"/45mm	Adolescent	TSG307R	TG307R	G46R	CG81
<b>MYA725LNA</b>	H205L-TH	7¼"	2¾"	L	453/1.0	1¾"/45mm	Female Adult Hand	TSG304L	TG304L	G45L	CG80
<b>MYA725RNA</b>	H206R-TH	7¼"	2¾"	R	453/1.0	1¾"/45mm	Female Adult Hand	TSG305R	TG305R	G46R	CG81
<b>MYA775LNA</b>	H207L-TH	7½"-7¾"	3"	L	493/1.1	2"/50mm	Female Adult Hand	TSG312L	TG312L	G47L	CG82
<b>MYA775RNA</b>	H208R-TH	7½"-7¾"	3"	R	493/1.1	2"/50mm	Female Adult Hand	TSG313R	TS313R	G48R	CG83
<b>MYA775LNA</b>	H207L-TH	7½"-7¾"	3"	L	493/1.1	2"/50mm	Male Adult Hand	TSG302L	TG302L	G47L	CG82
<b>MYA775RNA</b>	H208R-TH	7½"-7¾"	3"	R	493/1.1	2"/50mm	Male Adult Hand	TSG303R	TG303R	G48R	CG83
<b>MYA825LNA</b>	H209L-TH	7¾"-8¼"	3¼"	L	503/1.1	2⅛"/54mm	Male Adult Hand	TSG300L	TG300L	G49L	CG84
<b>MYA825RNA</b>	H210R-TH	7¾"-8¼"	3¼"	R	503/1.1	2⅛"/54mm	Male Adult Hand	TSG301R	TG301R	G50R	CG85

3/22/2010

\* Multi Control Plus™ Electric Hands

note: Silicone use "GSXXX"

**LTI** Liberating Technologies, Inc.

325 Hopping Brook Road, Suite A, Holliston, MA 01746-1456

www.liberatingtech.com

Phone 508-893-6363 FAX 508-893-9966