

LTI Data Sheets

Select Adult Foam Hands

FOAM HAND	hand	Side	Circumference	Gender	Length	Wrist dia.	Weight	Glove #	Glove #
part number	size (in)	L/R	mm w glove	finger-length	mm	mm	grams	silicone	PVC
FHL1	3½	left	212	male	140	54	228.6	TSGL1	TGL1
FHL2	3½	right	212	male	140	54	225.4	TSGL2	TGL2
FHL3	3¼	left	205	male	140	54	197.9	TSGL3	TGL3
FHL4	3¼	right	205	male	140	54	202.3	TSGL4	TGL4
FHL5	3	left	186	male/female	140	54	188.4	TSGL5	TGL5
FHL6	3	right	186	male/female	140	54	183.8	TSGL6	TGL6
FHL7	2¾	left	172	female	140	54	168.0	TSGL7	TGL7
FHL8	2¾	right	172	female	140	54	159.2	TSGL8	TGL8
FHL21	3¼	left	202	male-long	140	54	217.1	TSGL21	TGL9
FHL22	3¼	right	202	male-long	140	54	223.4	TSGL22	TGL10

Select glove color from RSLSteeper E-color swatches (E

Truefinish™ Glove lengths vary from 560-620 mm
FHL series Foam Hands require a (AP100) Wrist Stud

 Liberating Technologies, Inc.

325 Hopping Brook Road, Suite A, Holliston, MA 01746-1456

www.liberatingtech.com

Phone 508-893-6363 FAX 508-893-9966