

LTI Data Sheet Cosmetic Hands & Gloves

part #	RSL #	Metacarpal		side	description	Glove numbers			Wrist Options							
		width	circ.			LTI #	RSL #	#	Fixed round	Fixed oval	stud/plate	oval plate	friction*	OmniWrist	lightweight	connect pl

CHILDRENS HANDS

H301L	B11081	2"	5½"	L	Childs Foam Hand	G301L	19	40	WR150	NA	AP150	NA	MA163	VOMTH	WU150	WM150
H302R	B11082	2"	5½"	R	Childs Foam Hand	G302R	20	40	WR150	NA	AP150	NA	MA163	VOMTH	WU150	WM150
H303L	B11079	2¼"	6"	L	Childs Foam Hand	G303L	23	41	WR175	WO175	AP175	AP212	MA164	VOMTH	WU175	WM175
H304R	B11080	2¼"	6"	R	Childs Foam Hand	G304R	24	41	WR175	WO175	AP175	AP212	MA164	VOMTH	WU175	WM175
H305L	B11077	2¾"	6½"	L	Childs Foam Hand	G305L	17	42	WR175	WO175	AP175	AP212	MA164	VOMTH	WU175	WM175
H306R	B11078	2¾"	6½"	R	Childs Foam Hand	G306R	18	42	WR175	WO175	AP175	AP212	MA164	VOMTH	WU175	WM175
H307L	B11075	2½"	6¾"	L	Childs Foam Hand	G307L	15	43	WR175	WO175	AP175	AP212	MA165	VOMTH	WU175	WM175
H308R	B11076	2½"	6¾"	R	Childs Foam Hand	G308R	16	43	WR175	WO175	AP175	AP212	MA165	VOMTH	WU175	WM175
H309L	B11073	2½"	6¾"	L	Female Foam Hand	G309L	11	44	WR200	WO200	AP175	AP212	MA165	VOMTH	WU200	WE200
H310R	B11074	2½"	6¾"	R	Female Foam Hand	G310R	12	44	WR200	WO200	AP175	AP212	MA165	VOMTH	WU200	WE200

ADULT HANDS

New #	Old #	width	circ.	side	description	LTI #	RSL #		Fixed round	Fixed oval	stud/plate	oval plate	friction*	OmniWrist	lightweight	connect pl
FHL7	H315L	2¾"	7¼"	L	Female Foam Hand	TSGL7	L7	45	WR200	WO210	MA177	AP212	MA165	VOMTH	WU200	WE200
FHL8	H316R	2¾"	7¼"	R	Female Foam Hand	TSGL8	L8	45	WR200	WO210	MA177	AP212	MA165	VOMTH	WU200	WE200
FHL5	H317L	3"	7½"	L	Unisex Foam Hand	TSGL5	L5	48	WR210	WO210	MA177	AP212	MA165	VOMTH	WU210	WE210
FHL6	H318R	3"	7½"	R	Unisex Foam Hand	TSGL6	L6	48	WR210	WO210	MA177	AP212	MA165	VOMTH	WU210	WE210
FHL3	H321L	3¼"	7¾"	L	Mans Foam Hand	TSGL3	L3	50	WR210	WO210	MA177	AP212	MA165	VOMTH-05	WU210	WE210
FHL4	H322R	3¼"	7¾"	R	Mans Foam Hand	TSGL4	L4	50	WR210	WO210	MA177	AP212	MA165	VOMTH-05	WU210	WE210
FHL21	H323L	3¼"	7¾"	L	Mans long-finger Hand	TSGL21	L21	51	WR210	WO210	MA177	AP212	MA165	VOMTH-05	WU210	WE210
FHL22	H324R	3¼"	7¾"	R	Mans long-finger Hand	TSGL22	L22	51	WR210	WO210	MA177	AP212	MA165	VOMTH-05	WU210	WE210
FHL1	H325L	3½"	8"	L	Mans Foam Hand	TSGL1	L1	53	WR210	WO210	MA177	AP212	MA165	VOMTH-05	WU210	WE210
FHL2	H326R	3½"	8"	R	Mans Foam Hand	TSGL2	L2	53	WR210	WO210	MA177	AP212	MA165	VOMTH-05	WU210	WE210
H327L	<i>spec. order</i>	3¾"	8¼"	L	Mans Foam Hand	G27L	200	54	WR210	WO210	MA177	AP212	MA165	VOMTH-05	WU210	WE210
H328R	<i>spec. order</i>	3¾"	8¼"	R	Mans Foam Hand	G28R	201	54	WR210	WO210	MA177	AP212	MA165	VOMTH-05	WU210	WE210

* requires a wrist stud plate also

5/12/2010

LTI **Liberating Technologies, Inc.**

325 Hopping Brook Road, Suite A, Holliston, MA 01746-1456

www.liberatingtech.com

Phone 508-893-6363 FAX 508-893-9966