

Passive Prosthetic Hands for Adults and Children



- **Light-weight**
- **Cosmetically attractive**
- **Formable fingers/thumb**
- **Adult & children's sizes**
- **Gloves available in both men's and women's styles**
- **PVC, silicone and high-definition gloves available**

RSLSteeper Passive Foam Hands

RSLSteeper makes a series of passive foam hands for adults, teens, and children who require a purely cosmetic prosthetic hand. Intended for use with a cosmetic glove, these hands are not functional, but most can be shaped to the desired finger and thumb position. They consist of flexible foam on a formable plastic and metal skeleton, and hand plates are usually used to attach these foam hands to the socket. The hand dimensions specified indicate the width across the knuckles, in inches. The appropriate cosmetic glove is listed to the right of the product description. Always order at least one glove to accompany a passive hand, and specify the desired glove color. Order appropriate wrist from section 3.

For infants, we offer foam-filled passive hands from RSL Steeper and crawling hands from Centri. The foam-filled hands do *not* require a separate glove.

RSLSteeper Passive Hands for Adults and Children



H301L 2" Child's Foam Hand, left (shown)

H302R 2" Child's Foam Hand, right
measures 5.1" around palm
wrist-to-thumb 3"
wrist diameter 1.75"

glove G301L-SS

glove G302R-SS



H303L 2¼" Child's Foam Hand, left (shown)

H304R 2¼" Child's Foam Hand, right
measures 5.6" around palm
wrist-to-thumb 4"
wrist diameter 1.75"

glove G303L-SS

glove G304R-SS

Passive Prosthetic Hands for Adults and Children



H305L 2³/₈" Child's Foam Hand, left (shown)
H306R 2³/₈" Child's Foam Hand, right
measures 6.1" around palm
wrist-to-thumb 4"
wrist diameter 1.75"

glove G305L-SS
glove G306R-SS



H307L 2¹/₂" Child's Foam Hand, left (shown)
H308R 2¹/₂" Child's Foam Hand, right
measures 6.4" around palm
wrist-to-thumb 5"
wrist diameter 1.75"

glove G307L-SS
glove G308R-SS
Silicone
glove GS307L-SS
glove GS308R-SS



H309L 2¹/₂" Woman's Foam Hand, left (shown)
H310R 2¹/₂" Woman's Foam Hand, right
measures 6.6" around palm
wrist-to-thumb 5.5"
wrist diameter 2"

glove G309L-SS
glove G310R-SS
Silicone
glove GS309L-SS
glove GS310R-SS



H311L 2³/₄" Teen's Female Foam Hand, left
H312R 2³/₄" Teen's Female Foam Hand, right
H311L 2³/₄" Teen Foam Hand, left (shown)
H312R 2³/₄" Teen Foam Hand, right
measures 6.9" around palm
wrist-to-thumb 5.5"
wrist diameter 2.125"

glove G311L-SS
glove G312R-SS
glove G313L-SS
glove G314R-SS
Silicone
glove GS311L-SS
glove GS312R-SS
glove GS313L-SS
glove GS314R-SS



H315L 2³/₄" Woman's Foam Hand, left (shown)
H316R 2³/₄" Woman's Foam Hand, right
measures 7" around palm
wrist-to-thumb 5.5"
wrist diameter 2.125"

glove G315L-SS
glove G316R-SS
Silicone
glove GS315L-SS
glove GS316R-SS

Passive Prosthetic Hands for Adults and Children



H317L 3" Unisex Foam Hand, left (shown)

H318R 3" Unisex Foam Hand, right
measures 7.4" around palm
wrist-to-thumb 5.5"
wrist diameter 2.125"

glove G317L-SS
 glove G318R-SS
Silicone
 glove GS317L-SS
 glove GS318R-SS



H319L 3" Woman's Foam Hand, left (shown)

H320R 3" Woman's Foam Hand, right
measures 7.6" around palm
wrist-to-thumb 5.5"
wrist diameter 2.125"

glove G319L-SS
 glove G320R-SS
Silicone
 glove GS319L-SS
 glove GS320R-SS



H321L 3 1/4" Man's Foam Hand, short finger, left

H322R 3 1/4" Man's Foam Hand, short finger, right
measures 8.1" around palm
wrist-to-thumb 5.5"
wrist diameter 2.125"

glove G321L-SS
 glove G322R-SS
Silicone
 glove GS321L-SS
 glove GS322R-SS



H323L 3 1/4" Man's Foam Hand, long finger, left

H324R 3 1/4" Man's Foam Hand, long finger, right
measures 8.1" around palm
wrist-to-thumb 5.5"
wrist diameter 2.125"

glove G323L-SS
 glove G324R-SS
Silicone
 glove GS323L-SS
 glove GS324R-SS



H325L 3 1/2" Man's Foam Hand, left (shown)

H326R 3 1/2" Man's Foam Hand, right
measures 8.8" around palm
wrist-to-thumb 5.5"
wrist diameter 2.125"

glove G325L-SS
 glove G326R-SS
Silicone
 glove GS325L-SS
 glove GS326R-SS

Passive Prosthetic Hands for Adults and Children



H327L 3³/₄" Man's Foam Hand, left (shown)

H328R 3³/₄" Man's Foam Hand, right
measures 9.4" around palm
wrist-to-thumb 5.5"
wrist diameter 2.125"

glove G327L-SS
glove G328R-SS

Wrist options (see section 3 - WRISTS):

- NW** No Wrist
- RO** Round Wrist
- OV** Oval Wrist
- TH** Threaded Stud 1/2 -20
- LP** Lightweight Plate
- DP** Disconnect Plate
- LW** Lightweight



P80 Glove Lubricant, Pint
makes installing gloves easier, leaves no residue