

SKINERGY™ Leg Covers

SKINERGY™ protective leg covers are a new product manufactured using Tufsil 2, the latest innovation in silicone polymer technology. SKINERGY provides an extremely high level of cosmesis for the below knee or partial foot prosthesis, incorporating a unique seamless finish. Its improved tear resistance and 25% elasticity accommodates virtually all sockets and protects the prosthetic components of the leg. And each cover stretches to cover two foot sizes. SKINERGY is also stain resistant for longer life and can easily be cleaned with soap and water.

SKINERGY is now available in both split-toe and closed-toe versions making it suitable for all manufacturer's prosthetic feet and with two calf sizes, it can accommodate all users. Both styles are available in 12 colors and nails are hand painted for a "semi-custom" finish. Colors are selected from a color swatch sample. SKINERGY can be rolled down to access the underlying prosthetic hardware for maintenance and adjustments. To access the foot bolt for maintenance purposes, a 25 mm punch tool is available. And if the leg cover is damaged, a patch kit is offered consisting of color-matched patches and adhesive.

SKINERGY leg covers are available in four foot sizes and two calf sizes. The foot size is designated by the number (i.e. 22) which relates to the foot size in cm (22 cm or 220 mm). Each foot size has two calf sizes "small" or "large" designated by the trailing letter (S or L).



LC2XXX SKINERGY Protective Leg Cover, Split Toe
For use with all split-toe prosthetic feet

Part	Foot	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*	
Number	Size mm	minimum cm	maximum cm	minimum cm	maximum cm	minimum cm	maximum cm	Mid Calf cm	Cover Top cm
LC22LS	220/230	20	25	28	35	26	32.5	32	55
LC22RS	220/230	20	25	28	35	26	32.5	32	55
LC22LL	220/230	20	25	32	40	30	37.5	32	55
LC22RL	220/230	20	25	32	40	30	37.5	32	55
LC24LS	240/250	22	27.5	31	39	29	36	34.5	57
LC24RS	240/250	22	27.5	31	39	29	36	34.5	57
LC24LL	240/250	22	27.5	35	44	33	41	34.5	57
LC24RL	240/250	22	27.5	35	44	33	41	34.5	57
LC26LS	260/270	22	28	32	40	30	37.5	37	61
LC26RS	260/270	22	28	32	40	30	37.5	37	61
LC26LL	260/270	22	28	37	46	35	44	37	61
LC26RL	260/270	22	28	37	46	35	44	37	61
LC28LS	280/300	24	30	34	42.5	32	40	39.5	66
LC28RS	280/300	24	30	34	42.5	32	40	39.5	66
LC28LL	280/300	24	30	39	49	37	46	39.5	66
LC28RL	280/300	24	30	39	49	37	46	39.5	66

Closed-Toe SKINERGY Leg Covers



NS2XXX SKINERGY Protective Leg Cover, Closed Toe
For use with all closed-toe prosthetic feet

Part	Foot	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*	
Number	Size	minimum	maximum	minimum	maximum	minimum	maximum	Mid Calf	Cover Top
	mm	cm	cm	cm	cm	cm	cm	cm	cm
NS22LS	220/230	20	25	28	35	26	32.5	32	55
NS22RS	220/230	20	25	28	35	26	32.5	32	55
NS22LL	220/230	20	25	32	40	30	37.5	32	55
NS22RL	220/230	20	25	32	40	30	37.5	32	55
NS24LS	240/250	22	27.5	31	39	29	36	34.5	57
NS24RS	240/250	22	27.5	31	39	29	36	34.5	57
NS24LL	240/250	22	27.5	35	44	33	41	34.5	57
NS24RL	240/250	22	27.5	35	44	33	41	34.5	57
NS26LS	260/270	22	28	32	40	30	37.5	37	61
NS26RS	260/270	22	28	32	40	30	37.5	37	61
NS26LL	260/270	22	28	37	46	35	44	37	61
NS26RL	260/270	22	28	37	46	35	44	37	61
NS28LS	280/300	24	30	34	42.5	32	40	39.5	66
NS28RS	280/300	24	30	34	42.5	32	40	39.5	66
NS28LL	280/300	24	30	39	49	37	46	39.5	66
NS28RL	280/300	24	30	39	49	37	46	39.5	66

Dimensions are in cm, except foot size which is in mm

* Cover Length – from floor to mid-calf or cover top



P80 Cosmetic Cover Lubricant, Pint
makes installing leg covers easier, leaves no residue