

SKINERGY™ Leg Covers

ORDER FORM

SKINERGY split-toe (LC) and closed-toe (NS) leg covers are available in four foot sizes and two calf sizes. The Part Number designates the foot size which is stated in cm (22 cm = 220 mm). The side (L or R) is stated next. Each foot size has two calf sizes "small" and "large" designated by the trailing letter (S or L). Both styles are available in 12 colors. Specify Quantity of each size needed (left-hand column) and the colors (right-hand column) and FAX the form to LTI.

RSLSteeper Colors (see color swatch): 00, 0, 1, 3, 5, 7, 9, 11, 13, 15, 17, 19

Split-toe SKINERGY (LC#)

Quantity	Part	Foot	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*		Color #
	Number	Size mm	minimum cm	maximum cm	minimum cm	maximum cm	minimum cm	maximum cm	mid-calf cm	cover top cm	
	LC22LS	220/230	20	25	28	35	26	32.5	32	55	
	LC22RS	220/230	20	25	28	35	26	32.5	32	55	
	LC22LL	220/230	20	25	32	40	30	37.5	32	55	
	LC22RL	220/230	20	25	32	40	30	37.5	32	55	
	LC24LS	240/250	22	27.5	31	39	29	36	34.5	57	
	LC24RS	240/250	22	27.5	31	39	29	36	34.5	57	
	LC24LL	240/250	22	27.5	35	44	33	41	34.5	57	
	LC24RL	240/250	22	27.5	35	44	33	41	34.5	57	
	LC26LS	260/270	22	28	32	40	30	37.5	37	61	
	LC26RS	260/270	22	28	32	40	30	37.5	37	61	
	LC26LL	260/270	22	28	37	46	35	44	37	61	
	LC26RL	260/270	22	28	37	46	35	44	37	61	
	LC28LS	280/300	24	30	34	42.5	32	40	39.5	66	
	LC28RS	280/300	24	30	34	42.5	32	40	39.5	66	
	LC28LL	280/300	24	30	39	49	37	46	39.5	66	
	LC28RL	280/300	24	30	39	49	37	46	39.5	66	

Closed-toe SKINERGY (NS#)

Quantity	Part	Foot	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*		Color #
	Number	Size mm	minimum cm	maximum cm	minimum cm	maximum cm	minimum cm	maximum cm	mid-calf cm	cover top cm	
	NS22LS	220/230	20	25	28	35	26	32.5	32	55	
	NS22RS	220/230	20	25	28	35	26	32.5	32	55	
	NS22LL	220/230	20	25	32	40	30	37.5	32	55	
	NS22RL	220/230	20	25	32	40	30	37.5	32	55	
	NS24LS	240/250	22	27.5	31	39	29	36	34.5	57	
	NS24RS	240/250	22	27.5	31	39	29	36	34.5	57	
	NS24LL	240/250	22	27.5	35	44	33	41	34.5	57	
	NS24RL	240/250	22	27.5	35	44	33	41	34.5	57	
	NS26LS	260/270	22	28	32	40	30	37.5	37	61	
	NS26RS	260/270	22	28	32	40	30	37.5	37	61	
	NS26LL	260/270	22	28	37	46	35	44	37	61	
	NS26RL	260/270	22	28	37	46	35	44	37	61	
	NS28LS	280/300	24	30	34	42.5	32	40	39.5	66	
	NS28RS	280/300	24	30	34	42.5	32	40	39.5	66	
	NS28LL	280/300	24	30	39	49	37	46	39.5	66	
	NS28RL	280/300	24	30	39	49	37	46	39.5	66	

Dimensions are in cm, except foot size which is in mm

* Cover Length – from floor to mid-calf or cover top

Customer Name: _____ Date: _____

Customer Address: _____

Telephone Number: _____ P.O. # _____ Ship via: _____

UPS-ground is the default

FAX this form to Liberating Technologies 1-508-893-9966