

## SKINERGY PLUS™ Leg Covers

### ORDER FORM

SKINERGY PLUS split-toe (TLC) and closed-toe (TNS) leg covers are available in four foot sizes and two calf sizes. The Part Number designates the foot size which is stated in cm (22 cm = 220 mm). The side (L or R) is stated next. Each foot size has two calf sizes "small" and "large" designated by the trailing letter (S or L). Both styles are available in 19 colors. Specify Quantity of each size needed (left-hand column) and the colors (right-hand column) and FAX the form to LTI.

RSL Steeper Colors (see True Finish color swatches): E0,E1,E2,E3,E4,E5,E6,E7,E8,E9,E10,E11,E12,E13,E14,E15,E16,E17,E18

#### Split-toe SKINERGY PLUS(TLC#)

Quantity	Part	Foot	Ankle Circumference		Mid-Calf		Top of Cover		Cover Length*		Color #
	Number	Size mm	minimum cm	maximum cm	minimum cm	maximum cm	minimum cm	maximum cm	mid-calf cm	cover top cm	
	TLC22LS	220/230	20	25	28	35	26	32.5	32	55	
	TLC22RS	220/230	20	25	28	35	26	32.5	32	55	
	TLC22LL	220/230	20	25	32	40	30	37.5	32	55	
	TLC22RL	220/230	20	25	32	40	30	37.5	32	55	
	TLC24LS	240/250	22	27.5	31	39	29	36	34.5	57	
	TLC24RS	240/250	22	27.5	31	39	29	36	34.5	57	
	TLC24LL	240/250	22	27.5	35	44	33	41	34.5	57	
	TLC24RL	240/250	22	27.5	35	44	33	41	34.5	57	
	TLC26LS	260/270	22	28	32	40	30	37.5	37	61	
	TLC26RS	260/270	22	28	32	40	30	37.5	37	61	
	TLC26LL	260/270	22	28	37	46	35	44	37	61	
	TLC26RL	260/270	22	28	37	46	35	44	37	61	
	TLC28LS	280/300	24	30	34	42.5	32	40	39.5	66	
	TLC28RS	280/300	24	30	34	42.5	32	40	39.5	66	
	TLC28LL	280/300	24	30	39	49	37	46	39.5	66	
	TLC28RL	280/300	24	30	39	49	37	46	39.5	66	

#### Closed-toe SKINERGY PLUS(TNS#)

Quantity	Part	Foot	Ankle Circumference		Mid-Calf Circ.		Top of Cover Circ.		Cover Length*		Color #
	Number	Size mm	minimum cm	maximum cm	minimum cm	maximum cm	minimum cm	maximum cm	mid-calf cm	cover top cm	
	TNS22LS	220/230	20	25	28	35	26	32.5	32	55	
	TNS22RS	220/230	20	25	28	35	26	32.5	32	55	
	TNS22LL	220/230	20	25	32	40	30	37.5	32	55	
	TNS22RL	220/230	20	25	32	40	30	37.5	32	55	
	TNS24LS	240/250	22	27.5	31	39	29	36	34.5	57	
	TNS24RS	240/250	22	27.5	31	39	29	36	34.5	57	
	TNS24LL	240/250	22	27.5	35	44	33	41	34.5	57	
	TNS24RL	240/250	22	27.5	35	44	33	41	34.5	57	
	TNS26LS	260/270	22	28	32	40	30	37.5	37	61	
	TNS26RS	260/270	22	28	32	40	30	37.5	37	61	
	TNS26LL	260/270	22	28	37	46	35	44	37	61	
	TNS26RL	260/270	22	28	37	46	35	44	37	61	
	TNS28LS	280/300	24	30	34	42.5	32	40	39.5	66	
	TNS28RS	280/300	24	30	34	42.5	32	40	39.5	66	
	TNS28LL	280/300	24	30	39	49	37	46	39.5	66	
	TNS28RL	280/300	24	30	39	49	37	46	39.5	66	

Dimensions are in cm, except foot size which is in mm

\* Cover Length – from floor to mid-calf or cover top

Customer Name: \_\_\_\_\_ Date: \_\_\_\_\_

Customer Address: \_\_\_\_\_

Telephone Number: \_\_\_\_\_ P.O. # \_\_\_\_\_ Ship via: \_\_\_\_\_

UPS-ground is the default

FAX this form to Liberating Technologies 1-508-893-9966