

**SKINERGY™ Leg Covers**  
Conversion Table

**Split-toe SKINERGY (LC#)**

Part	Foot	Foot	Ankle Circumference				Mid-Calf Circumference				Top of Cover Circumference				Cover Length*			
Number	Size	Size	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	Mid Calf	Mid Calf in	Cover Top	Cover Top
#	mm	inches	cm	cm	in	in	cm	cm	in	in	cm	cm	in	in	cm	in	cm	in
LC22LS	220/230	8.7/9.1	20	25	7.9	9.8	28	35	11.0	13.8	26	32.5	10.2	12.8	32	12.6	55	21.7
LC22RS	220/230	8.7/9.1	20	25	7.9	9.8	28	35	11.0	13.8	26	32.5	10.2	12.8	32	12.6	55	21.7
LC22LL	220/230	8.7/9.1	20	25	7.9	9.8	32	40	12.6	15.7	30	37.5	11.8	14.8	32	12.6	55	21.7
LC22RL	220/230	8.7/9.1	20	25	7.9	9.8	32	40	12.6	15.7	30	37.5	11.8	14.8	32	12.6	55	21.7
LC24LS	240/250	9.4/9.8	22	27.5	8.7	10.8	31	39	12.2	15.4	29	36	11.4	14.2	34.5	13.6	57	22.4
LC24RS	240/250	9.4/9.8	22	27.5	8.7	10.8	31	39	12.2	15.4	29	36	11.4	14.2	34.5	13.6	57	22.4
LC24LL	240/250	9.4/9.8	22	27.5	8.7	10.8	35	44	13.8	17.3	33	41	13.0	16.1	34.5	13.6	57	22.4
LC24RL	240/250	9.4/9.8	22	27.5	8.7	10.8	35	44	13.8	17.3	33	41	13.0	16.1	34.5	13.6	57	22.4
LC26LS	260/270	10.2/10.6	22	28	8.7	11.0	32	40	12.6	15.7	30	37.5	11.8	14.8	37	14.6	61	24.0
LC26RS	260/270	10.2/10.6	22	28	8.7	11.0	32	40	12.6	15.7	30	37.5	11.8	14.8	37	14.6	61	24.0
LC26LL	260/270	10.2/10.6	22	28	8.7	11.0	37	46	14.6	18.1	35	44	13.8	17.3	37	14.6	61	24.0
LC26RL	260/270	10.2/10.6	22	28	8.7	11.0	37	46	14.6	18.1	35	44	13.8	17.3	37	14.6	61	24.0
LC28LS	280/300	11.0/11.8	24	30	9.4	11.8	34	42.5	13.4	16.7	32	40	12.6	15.7	39.5	15.6	66	26.0
LC28RS	280/300	11.0/11.8	24	30	9.4	11.8	34	42.5	13.4	16.7	32	40	12.6	15.7	39.5	15.6	66	26.0
LC28LL	280/300	11.0/11.8	24	30	9.4	11.8	39	49	15.4	19.3	37	46	14.6	18.1	39.5	15.6	66	26.0
LC28RL	280/300	11.0/11.8	24	30	9.4	11.8	39	49	15.4	19.3	37	46	14.6	18.1	39.5	15.6	66	26.0

**Closed-toe SKINERGY (NS#)**

Part	Foot Size	Foot Size	Ankle Circumference				Mid-Calf Circumference				Top of Cover Circumference				Cover Length*			
Number	mm	inches	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	Mid Calf	Mid Calf in	Cover Top	Cover Top
#	mm	inches	cm	cm	in	in	cm	cm	in	in	cm	cm	in	in	cm	in	cm	in
NS22LS	220/230	8.7/9.1	20	25	7.9	9.8	28	35	11.0	13.8	26	32.5	10.2	12.8	32	12.6	55	21.7
NS22RS	220/230	8.7/9.1	20	25	7.9	9.8	28	35	11.0	13.8	26	32.5	10.2	12.8	32	12.6	55	21.7
NS22LL	220/230	8.7/9.1	20	25	7.9	9.8	32	40	12.6	15.7	30	37.5	11.8	14.8	32	12.6	55	21.7
NS22RL	220/230	8.7/9.1	20	25	7.9	9.8	32	40	12.6	15.7	30	37.5	11.8	14.8	32	12.6	55	21.7
NS24LS	240/250	9.4/9.8	22	27.5	8.7	10.8	31	39	12.2	15.4	29	36	11.4	14.2	34.5	13.6	57	22.4
NS24RS	240/250	9.4/9.8	22	27.5	8.7	10.8	31	39	12.2	15.4	29	36	11.4	14.2	34.5	13.6	57	22.4
NS24LL	240/250	9.4/9.8	22	27.5	8.7	10.8	35	44	13.8	17.3	33	41	13.0	16.1	34.5	13.6	57	22.4
NS24RL	240/250	9.4/9.8	22	27.5	8.7	10.8	35	44	13.8	17.3	33	41	13.0	16.1	34.5	13.6	57	22.4
NS26LS	260/270	10.2/10.6	22	28	8.7	11.0	32	40	12.6	15.7	30	37.5	11.8	14.8	37	14.6	61	24.0
NS26RS	260/270	10.2/10.6	22	28	8.7	11.0	32	40	12.6	15.7	30	37.5	11.8	14.8	37	14.6	61	24.0
NS26LL	260/270	10.2/10.6	22	28	8.7	11.0	37	46	14.6	18.1	35	44	13.8	17.3	37	14.6	61	24.0
NS26RL	260/270	10.2/10.6	22	28	8.7	11.0	37	46	14.6	18.1	35	44	13.8	17.3	37	14.6	61	24.0
NS28LS	280/300	11.0/11.8	24	30	9.4	11.8	34	42.5	13.4	16.7	32	40	12.6	15.7	39.5	15.6	66	26.0
NS28RS	280/300	11.0/11.8	24	30	9.4	11.8	34	42.5	13.4	16.7	32	40	12.6	15.7	39.5	15.6	66	26.0
NS28LL	280/300	11.0/11.8	24	30	9.4	11.8	39	49	15.4	19.3	37	46	14.6	18.1	39.5	15.6	66	26.0
NS28RL	280/300	11.0/11.8	24	30	9.4	11.8	39	49	15.4	19.3	37	46	14.6	18.1	39.5	15.6	66	26.0