

RSL Steeper Mechanical Hands



VC Dorsal pull adult and child's hand

- **Voluntary Close (VC)**
- **Two sizes for children, four sizes for adults**
- **Positive lock on larger hands for carrying objects**
- **Adjustable grip force**
- **Dorsal and palmar pull cords**
- **Cosmetic gloves in 11 colors**

These mechanical hands are designed for control with either a palmar (medial) or a dorsal (lateral) pull cord, ensuring that they can be fitted to most prosthetic systems. They are also available in both Voluntary Open (VO) versions. Order appropriate wrist from section 3.

Voluntary Close Hands



H57VCL 2" Dorsal Pull Hand VC, left
H58VCR 2" Dorsal Pull Hand VC, right (shown)
measures 5" around palm, 1½" diameter wrist
use with AP150 Hand Connector Plate, see section 3

glove G29L-SS
 glove G30R-SS



H61VCL 2¼" Dorsal Pull Hand VC, left (shown)
H62VCR 2¼" Dorsal Pull Hand VC, right
measures 5¾" around palm, 1¾" diameter wrist
use with AP175 Hand Connector Plate, see section 3

glove G31L-SS
 glove G32R-SS



H65VCL 2½" Dorsal Pull Hand VC, left (shown)
H66VCR 2½" Dorsal Pull Hand VC, right
measures 6¼" around palm, 1¾" diameter wrist
use with AP175 Hand Connector Plate, see section 3

glove G33L-SS
 glove G34R-SS

RSLSteeper Mechanical Hands



H69VCL 2³/₄" Dorsal Pull Hand VC, left (shown)
H70VCR 2³/₄" Dorsal Pull Hand VC, right
measures 7" around palm, 2" diameter wrist
use with AP200 Hand Connector Plate, see section 3

glove G35L-SS
 glove G36R-SS



H73VCL 3" Dorsal Pull Hand VC, left (shown)
H74VCR 3" Dorsal Pull Hand VC, right
measures 7" around palm, 2" diameter wrist
use with AP200 Hand Connector Plate, see section 3

glove G37L-SS
 glove G38R-SS



H77VCL 3¹/₄" Dorsal Pull Hand VC, left
H78VCR 3¹/₄" Dorsal Pull Hand VC, right (shown)
measures 7¹/₂" around palm, 2" diameter wrist
use with AP200 Hand Connector Plate, see section 3

glove G39L-SS
 glove G40R-SS

Cosmetic Glove Colors:



When ordering RSLSteeper gloves, replace the **-SS** suffix with a color code from the list. Some uncommon colors may need to be special ordered, taking 3-4 weeks.

For proper color match, request color swatches or use Otto Bock colors.

Unless otherwise stated Colour 3 will be supplied.