

## ***RSL Steeper Mechanical Hands & Hooks***



*VO Dorsal and Palmar pull adult and child's hands*

- **Voluntary Open (VO)**
- **Two sizes for children, four sizes for adults**
- **Positive lock on larger sizes for carrying objects**
- **Adjustable grip force**
- **Dorsal and palmar pull cords**
- **Cosmetic gloves in 11 colors**

These Voluntary Open mechanical hands are designed for control with either a palmar (medial) or a dorsal (lateral) pull cord, ensuring that they can be fitted to most prosthetic systems. They are also available in both Voluntary Close (VC) versions. Order appropriate wrist from section 3.

### ***Voluntary Open Hands***



**H57L** 2" Palmar Pull Hand VO, left (shown)  
**H58R** 2" Palmar Pull Hand VO, right  
*measures 5" around palm, 1½" diameter wrist*  
*use with AP150 Hand Connector Plate, see section 3*

glove G29L-SS  
 glove G30R-SS



**H59L** 2" Dorsal Pull Hand VO, left (shown)  
**H60R** 2" Dorsal Pull Hand VO, right  
*measures 5" around palm, 1½" diameter wrist*  
*use with AP150 Hand Connector Plate, see section 3*

glove G29L-SS  
 glove G30R-SS



**H61L** 2¼" Palmar Pull Hand VO, left (shown)  
**H62R** 2¼" Palmar Pull Hand VO, right  
*measures 5¾" around palm, 1¾" diameter wrist*  
*use with AP175 Hand Connector Plate, see section 3*

glove G31L-SS  
 glove G32R-SS

## ***RSL Steeper Mechanical Hands & Hooks***



**H63L** 2¼" Dorsal Pull Hand VO, left (shown)  
**H64R** 2¼" Dorsal Pull Hand VO, right  
*measures 5¾" around palm, 1¾" diameter wrist*  
*use with AP175 Hand Connector Plate, see section 3*

glove G31L-SS  
glove G32R-SS



**H65L** 2½" Palmar Pull Hand VO, left (shown)  
**H66R** 2½" Palmar Pull Hand VO, right  
*measures 6¼" around palm, 1¾" diameter wrist*  
*use with AP175 Hand Connector Plate, see section 3*

glove G33L-SS  
glove G34R-SS



**H67L** 2½" Dorsal Pull Hand VO, left (shown)  
**H68R** 2½" Dorsal Pull Hand VO, right  
*measures 6¼" around palm, 1¾" diameter wrist*  
*use with AP175 Hand Connector Plate, see section 3*

glove G33L-SS  
glove G34R-SS



**H69L** 2¾" Palmar Pull Hand VO, left (shown)  
**H70R** 2¾" Palmar Pull Hand VO, right  
*measures 7" around palm, 2" diameter wrist*  
*use with AP200 Hand Connector Plate, see section 3*

glove G35L-SS  
glove G36R-SS



**H71L** 2¾" Dorsal Pull Hand VO, left (shown)  
**H72R** 2¾" Dorsal Pull Hand VO, right  
*measures 7" around palm, 2" diameter wrist*  
*use with AP200 Hand Connector Plate, see section 3*

glove G35L-SS  
glove G36R-SS

## ***RSLSteeper Mechanical Hands & Hooks***



**H73L** 3" Palmar Pull Hand VO, left (shown)  
**H74R** 3" Palmar Pull Hand VO, right  
*measures 7" around palm, 2" diameter wrist*  
*use with AP200 Hand Connector Plate, see section 3*

glove G37L-SS  
glove G38R-SS



**H75L** 3" Dorsal Pull Hand VO, left (shown)  
**H76R** 3" Dorsal Pull Hand VO, right  
*measures 7" around palm, 2" diameter wrist*  
*use with AP200 Hand Connector Plate, see section 3*

glove G37L-SS  
glove G38R-SS



**H77L** 3/4" Palmar Pull Hand VO, left (shown)  
**H78R** 3/4" Palmar Pull Hand VO, right  
*measures 7 1/2" around palm, 2" diameter wrist*  
*use with AP200 Hand Connector Plate, see section 3*

glove G39L-SS  
glove G40R-SS



**H79L** 3/4" Dorsal Pull Hand VO, left (shown)  
**H80R** 3/4" Dorsal Pull Hand VO, right  
*measures 7 1/2" around palm, 2" diameter wrist*  
*use with AP200 Hand Connector Plate, see section 3*

glove G39L-SS  
glove G40R-SS

Wrist options (see Section 3 - WRISTS):

- NW** No Wrist
- RO** Round Wrist
- OV** Oval Wrist
- TH** Threaded Stud 1/2 -20
- LP** Lightweight Plate
- DP** Disconnect Plate
- LW** Lightweight

# *RSL Steeper Mechanical Hands & Hooks*

## *Light-Weight Split-Hooks*



**TH201L** Titanium Split Hook, Left  
Weighs only 132 grams, 1/2-20 UNF thread



**TH202R** Titanium Split Hook, Right  
Weighs only 132 grams, 1/2-20 UNF thread