

Prosthetic Wrists

Wrists for VASI Electric Hands

The VASI Short Oval and Short Round wrists are suitable for children with long residual limbs and natural pronation. The Short Friction wrist is suitable for children with a long limb needing a friction adjustment. The standard Friction wrist and the LTI-Houston FlexiWrist are suitable for shorter limbs.



VSO03 Short Oval Wrist 0-3 for VASI 0-3 Electric Hands (shown)

VSO26 Short Oval Wrist 2-6 for VASI 2-6 Electric Hands

VSO59 Short Oval Wrist 5-9 for VASI 5-9 & 7-11 Electric Hands



VSR03 Short Round Wrist 0-3 for VASI 0-3 Electric Hands (shown)

Single groove for lamination; 1.38" (35mm) in diameter

VSR26 Short Round Wrist 2-6 for VASI 2-6 Electric Hands

Single groove for lamination; 1.46" (37mm) in diameter

VSR59 Short Round Wrist 5-9 for VASI 5-9 & 7-11 Electric Hands

Single groove for lamination; 1.65" (42mm) in diameter



VSF03 Short Friction Wrist 0-3 for VASI 0-3 Electric Hands (shown)

Single groove for lamination; 1.46" (37mm) in diameter

VSF26 Short Friction Wrist 2-6 for VASI 2-6 Electric Hands

Single groove for lamination; 1.70" (43mm) in diameter

VSF59 Short Friction Wrist 5-9 for VASI 5-9 & 7-11 Electric Hands

Single groove for lamination; 1.70" (43mm) in diameter



VFR03 Friction Wrist 0-3 for VASI 0-3 Electric Hands (shown)

Three-groove lamination collar, 1.13" (28.6mm) in diameter

VFR26 Friction Wrist 2-6 for VASI 2-6 Electric Hands

Three-groove lamination collar, 1.46" (37mm) in diameter

VFR59 Friction Wrist 5-9 for VASI 5-9 & 7-11 Electric Hands

Three-groove lamination collar, 1.46" (37mm) in diameter



VEF03 Elbow Friction Wrist 0-3 for VASI 0-3 Electric Hands

VEF26 Elbow Friction Wrist 2-6 for VASI 2-6 Electric Hands

VEF59 Elbow Friction Wrist 5-9 for VASI 5-9 & 7-11 Electric Hands

Wide three-groove collar for cementing into VASI elbow forearm 1.57" (40mm) in diameter