

SKINERGY™

The award winning
lower limb silicone cosmesis

Split toe



* INDEPENDENT LIVING DESIGN
AWARD FOR INNOVATION

SKINERGY™ - The award

Following international success of the award winning SKINERGY cosmesis for **split toe** foot designs, **RSLSTEEPER** is delighted to announce the arrival of SKINERGY for **closed toe** foot types.

Manufactured using Tufsil 2, the latest innovation in silicone polymer technology, SKINERGY provides an extremely high level of cosmesis for the below knee user, incorporating a unique seamless finish. Furthermore the product provides a level of protection for the hardware contained within the prosthesis.

Avoiding the tacky feel of other silicone products, the unique properties of SKINERGY include significantly improved tear resistance together with up to 25% stretch, allowing even bulbous sockets to be accommodated. These characteristics provide additional benefits, with each SKINERGY comfortably stretching to fit two foot sizes and variations in foot types.

SKINERGY is simply rolled on and off to allow access to underlying hardware for maintenance and adjustment, and can also be easily cleaned using soap and water.

In response to clinical feedback, we have also lengthened both ranges by approximately 20% to provide an even smoother transition over the proximal socket brim. SKINERGY for both **split toe** and **closed toe** designs is available in a range of 12 color shades, and all nails are hand painted to complete the almost custom finish.

To access the foot bolt for maintenance purposes if required, a punch tool (NCT9740) is available to create an access hole. The Color Patch Kit (E22630) includes adhesive and color matched silicone patches to provide a suitable seal. Spare patches are also available (E22632).

SKINERGY™ Cosmesis benefits include:

- **Advanced silicone material - durable and stain resistant**
- **Assists in protection of the prosthesis**
- **Easily removable for hardware adjustment, maintenance and repair**
- **Available for any foot with a split or closed toe design - sizes 220mm to 300mm**
- **Off the shelf availability**
- **Easy to clean using soap and water**
- **Custom nail finish**
- **Unique seamless design**
- **Available in 12 skin tones**
- **A true silicone cosmesis, not just a cosmetic cover**



skinergy
skinergy
skinergy
skinergy
skinergy
skinergy

ard winning lower lim

ORDERING INFORMATION

Part No Split Toe	Part No Closed Toe	Foot Size mm	Ankle Circ.		Mid Calf Circ.		Top of Cover Circ.		Cover Length from Floor to - Mid Calf Cover Top cm	
			Min cm	Max cm	Min cm	Max cm	Min cm	Max cm		
LC22LS	NS22LS	220 / 230	20	25	28	35	26	32.5	32	55
LC22RS	NS22RS	220 / 230	20	25	28	35	26	32.5	32	55
LC22LL	NS22LL	220 / 230	20	25	32	40	30	37.5	32	55
LC22RL	NS22RL	220 / 230	20	25	32	40	30	37.5	32	55
LC24LS	NS24LS	240 / 250	22	27.5	31	39	29	36	34.5	57
LC24RS	NS24RS	240 / 250	22	27.5	31	39	29	36	34.5	57
LC24LL	NS24LL	240 / 250	22	27.5	35	44	33	41	34.5	57
LC24RL	NS24RL	240 / 250	22	27.5	35	44	33	41	34.5	57
LC26LS	NS26LS	260 / 270	22	28	32	40	30	37.5	37	61
LC26RS	NS26RS	260 / 270	22	28	32	40	30	37.5	37	61
LC26LL	NS26LL	260 / 270	22	28	37	46	35	44	37	61
LC26RL	NS26RL	260 / 270	22	28	37	46	35	44	37	61
LC28LS	NS28LS	280 / 300	24	30	34	42.5	32	40	39.5	66
LC28RS	NS28RS	280 / 300	24	30	34	42.5	32	40	39.5	66
LC28LL	NS28LL	280 / 300	24	30	39	49	37	46	39.5	66
LC28RL	NS28RL	280 / 300	24	30	39	49	37	46	39.5	66

DESCRIPTION

Skinergy™ Color Patch Kit

State color required e.g. E22630/3

(Includes adhesives & 5 patches)

Part No

E22630

Skinergy™ Color Swatch

E18733

Skinergy™ 25mm Punch Tool

NCT9740

Skinergy™ Patches (5)

State color required e.g. E22632/3

E22632



skinenergy silicone cosmesis

skinenergy
skinenergy
skinenergy
skinenergy
skinenergy
skinenergy
skinenergy
skinenergy
skinenergy

Closed toe



www.skinenergy.info